

Plus = Maximum effectiveness: The new acupuncture needle Aktiv Plus

A small tool with a lot of potential

The specific needles used in a treatment are just as important to success as the skill of the acupuncturist, and the demands placed on needles are high. Good acupuncture needles should penetrate a patient's skin causing as little discomfort as possible. In order to achieve this, great emphasis is given to the design of a needle's tip and the smoothness of its shaft. The problem, however, is that smooth needle shafts lessen the desired effect, particularly in areas that would benefit from a high level of stimulation. Because of this, Schwa-Medico – a company with decades of experience in developing products for and working together with acupuncturists – set out to place particular emphasis on the needle shaft from a medical-technology perspective and to develop an acupuncture needle that would offer a higher level of stimulation and a tip that would be finished to the highest specifications. The needle glides smoothly through the skin and provides a high level of stimulation at the relevant points. All of these specifications were successfully met by developing a needle with a special shaft with spindle-shaped grooves that meet at a perfectly smooth tip. The new acupuncture needle Singer Aktiv Plus therefore offers totally new possibilities in acupuncture treatments.

A needle with a high profile


High quality steel, a perfectly formed tip, electropolished three times and a highly proven ergonomic handle round off the profile of the Singer Aktiv Plus. The high-specification tip reduces discomfort when penetrating the skin and makes skin penetration significantly easier. The spindle-shaped shaft increases the needle's stimulatory effect and ensures safe and patient-friendly acupuncture. An additional benefit: It also helps to trigger the sensation of De-Qi in the outer layers of skin and the needle's handle offers maximum comfort when being handled.



At one glance

Characteristics	Advantages
<ul style="list-style-type: none">■ spindle shaped needle shaft	<ul style="list-style-type: none">■ increases the needle's stimulatory effect■ ensures safe and patient-friendly acupuncture■ helps to trigger the sensation of De-Qi in the outer layers of skin
<ul style="list-style-type: none">■ perfectly formed tip	<ul style="list-style-type: none">■ reduces discomfort when penetrating the skin■ makes skin penetration significantly easier
<ul style="list-style-type: none">■ ergonomic needle handle	<ul style="list-style-type: none">■ offers maximum comfort when being handled

In comparison



Clearly visible in enlargement

The spindle shaped needle shaft of the Singer Aktiv plus

Optimum stimulation -

The new acupuncture needle
by schwa-medico



schwa-medico
Medizinische Apparate
Vertriebsgesellschaft mbH
Export Department
Gehrstraße 4
D-35630 Ehringshausen
Tel. + 49 6443 8333-717
Fax + 49 6443 8333-766
E-Mail friedrich@schwa-medico.de
www.schwa-medico.com

